



596PR-01

G-Biosciences ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ [technical@GBiosciences.com](mailto:technical@GBiosciences.com)

A Geno Technology, Inc. (USA) brand name

# Silica Magnetic Beads

(Cat. # 786-915, 786-916, 786-917)



think proteins! think G-Biosciences [www.GBiosciences.com](http://www.GBiosciences.com)

INTRODUCTION ..... 3

ITEMS SUPPLIED ..... 3

STORAGE CONDITIONS ..... 3

SPECIFICATIONS ..... 3

PRECAUTIONS ..... 3

PROTOCOL ..... 4

    ADDITIONAL ITEMS REQUIRED ..... 4

    PREPARATION OF SILICA MAGNETIC BEADS ..... 4

    PURIFICATION OF NUCLEIC ACID ..... 4

RELATED PRODUCTS ..... 5

## INTRODUCTION

G-Biosciences Silica Magnetic Beads are  $\text{Fe}_3\text{O}_4$  magnetic beads coated with a silicon dioxide ( $\text{SiO}_2$ ) layer. Since silica is able to bind to the nucleic acids, G-Biosciences Silica Magnetic Beads serve as a simple and efficient tool for plasmid DNA purification for transfection or sequencing applications, genomic DNA purification for research or clinical applications, RNA purification for qPCR analysis, or PCR product clean-up for downstream analysis.

## ITEMS SUPPLIED

Cat. #	Description	Size
786-915	Silica Magnetic Beads	5ml
786-916	Silica Magnetic Beads	25ml
786-917	Silica Magnetic Beads	100ml

## STORAGE CONDITIONS

The beads are shipped at ambient temperature. Upon arrival, store the beads at 4°C. If stored and handled correctly the beads have a 1 year shelf life.

## SPECIFICATIONS

$\text{Fe}_3\text{O}_4$  beads coated with silicon dioxide ( $\text{SiO}_2$ ) of an average 2.5-4.5 $\mu\text{m}$  in diameter for the binding of nucleic acids. Binding capacity is 4mg DNA/ml beads. G-Biosciences Silica Magnetic Beads are supplied in phosphate buffered saline, pH 7.4 with 0.09% Sodium Azide and 0.02% Tween-20.

## PRECAUTIONS

- Do not freeze the magnetic beads
- Do not store near magnetic sources

## PROTOCOL

### ***Additional Items Required***

- Binding Buffer: 4M guanidine thiocyanate, 40mM Tris, 17.6mM EDTA, pH 8.0
- Wash Buffer: 10mM Tris-HCl, 1mM EDTA, 70% ethanol, pH8.0
- Elution Buffer: TE Buffer (10mM Tris-HCl, 1mM EDTA, pH8.0)
- Magnetic Stand or magnet

### ***Preparation of Silica Magnetic Beads***


1. Resuspend G-Biosciences Silica Magnetic Beads thoroughly by pipetting or vortex the vial.
2. Transfer adequate amount of beads into a clean tube.
3. Place the tube on the magnetic stand for 30-60 seconds to immobilize the beads at tube wall.
4. Discard the supernatant by aspiration with a pipette.
5. Remove the tube from magnetic stand.
6. Add 100  $\mu$ l Elution Buffer (or ddH<sub>2</sub>O) and resuspend the beads by pipetting or vortex.
7. Place the tube on the magnetic stand for 30-60 seconds to immobilize the beads at tube wall.
8. Discard the supernatant, and then remove the tube from the magnetic stand.
9. Repeat steps 6-8 twice.

### ***Purification of Nucleic Acid***

1. Mix 10 $\mu$ l sample and 90 $\mu$ l Binding Buffer with magnetic beads thoroughly by pipetting.
2. Incubate with tilt rotation for 2 minutes at room temperature.
3. Place the tube on the magnetic stand for 30-60 seconds to immobilize the beads at tube wall.
4. Discard (or collect) the supernatant as unbound substances by aspiration with a pipette, and then remove the tube from the magnetic stand.
5. Add 100 $\mu$ l Wash Buffer and resuspend the beads by pipetting.
6. Place the tube on the magnetic stand for 30-60 seconds to immobilize the beads at tube wall.
7. Discard (or collect) the supernatant as unbound substances, and then remove the tube from the magnetic stand.
8. Repeat steps 5-7 twice.
9. Air-dry for 5-20 min.
10. Add 10-100 $\mu$ l Elution Buffer (or ddH<sub>2</sub>O) and resuspend the beads complex by vortex or shaking.
11. Incubate with tilt rotation for 3 minutes at room temperature.
12. Place the tube on the magnetic stand for 30-60 seconds and collect the supernatant to a clean tube.

## RELATED PRODUCTS

Download our Molecular Biology Handbook.



<http://info.gbiosciences.com/complete-molecular-biology-handbook>

For other related products, visit our website at [www.GBiosciences.com](http://www.GBiosciences.com) or contact us.

Last saved: 5/19/2015 CMH

*This page is intentionally left blank*

*This page is intentionally left blank*



[www.GBiosciences.com](http://www.GBiosciences.com)