

## ProCipitate™

### *Superior Substitute to Phenol/Chloroform for DNA & RNA Isolation*



- Removes only the contaminants & leaves DNA alone
- Improves yield of DNA over alternative bind and elute systems
- Adaptable to any sample size, and can be automated
- Key component of the ProPrep™ line of application specific kits

ProCipitate™ is a unique protein extraction reagent based upon patented elastomeric polyelectrolytes. The polymer chains of ProCipitate™ are initially extended in a high energy state due to an overall net negative charge. When introduced to protein solutions, the charges are neutralized and the polymer chains collapse to a more favorable energy state; DNA and RNA remain unreacted.

ProCipitate™ is non-hazardous and can replace phenol/chloroform with the additional benefits of solid-phase suspensions: adaptability to filtration and automation. It is routinely used for Plasmids, Cosmids, BACs, and Genomic DNA, as well as RNA. ProCipitate™ can also be used to remove Proteinase K and other enzymes.

ProCipitate™ provides high quality DNA suitable for automated sequencing, Southern blotting, and restriction digestion. ProCipitate™ is available as a suspension reagent and in ProPrep™ kits for specific applications and high-throughput 96 well filter formats.

Sample Size	ProCipitate™ Typical Usage
10 ml Yeast Culture	1-2 ml
Mouse Tail	250-500 µl
4 mm Plant Leaf	100-200 µl





Blanket order or large quantity discounts are available for any of these items.

Product Description	Size & Quantity	Item Number	Price
ProCipitate™	30 ml	P0050-30	\$335
ProCipitate™	100 ml	P0050-100	\$870
ProCipitate™	50 ml	P0050-500	\$3206.67
<b>High-Throughput Genomic DNA Isolation Kits</b>			
ProPrep™ Genomic 96	96, 50µl whole blood	PPG-96	\$280
ProPrep™ Genomic 960	10x96, 50µl whole blood	PPG-960	\$1980
ProPrep™ Genomic SM-50	50, 0.1ml whole blood	PBKSM-50	\$210
ProPrep™ Genomic XL-2	2, 10ml whole blood	PBKXL-2	\$235
ProPrep™ Genomic XL-10	10, 10ml whole blood	PBKXL-10	\$665
ProPrep™ Genomic Blood Card 96	96 blood spots	PBC-96	\$180
ProPrep™ Genomic Blood Card 960	10x96 blood spots	PBC-960	\$1441

Product Description	Size & Quantity	DNA Yield* Per Prep	Item Number	Price
<b>High-Throughput BAC &amp; Plasmid DNA Isolation</b>				
ProPrep™ BAC Mini 100	100, 2 ml cultures	2 µg BAC	PMK-100	\$255
ProPrep™ BAC Mini 1000	1000, 2 ml cultures	2 µg BAC	PMK-1000	\$1078
ProPrep™ BAC 96	1x96, 2 ml cultures	2 µg BAC	PLF-96	\$310
ProPrep™ BAC 960	10x96, 2 ml cultures	2 µg BAC	PLF-960	1595
ProPrep™ BAC Omni 200	Up to 200 ml cultures	50 µg BAC	PPO-200	\$450
ProPrep™ BAC Omni 1000	Up to 1000 ml cultures	250 µg BAC	PPO-1000	\$1078
ProPrep™ Plasmid 4x96 (384 preps)	4x96, 250µl cultures	2-3 µg, ≈ 150 ng/µl	PPF-4x96	\$385
ProPrep™ Plasmid 60x96 (5760 preps)	60x96, 250µl cultures	2-3 µg, ≈ 150 ng/µl	PPF-60x96	\$4800

\*Typical under recommended growth conditions, actual results may vary.

## STORAGE

ProCipitate™ is an aqueous suspension polyelectrolyte in distilled water. Shake well before use. The reagent when not used must be kept sealed and stored at 4°C. ProCipitate™ retains full activity when stored at 4°C for about 6 months. For long term storage, please contact technical services. ProCipitate™ performs optimally in a pH range of approximately 4 to 6, however the polyelectrolyte is sufficiently acidic (pH 4) to lower the pH to within its optimal working range in most applications.



BIOTECH SUPPORT GROUP

## Performance Characteristics Typical Volume Usage

Protein	ProCipitate™: Sample	Removal
BSA, PBS @ 30 mg/ml	1 : 1	>99%
BSA, 1%SDS @ 30 mg/ml	1 : 1	>99%
BSA, 3M GuSCN @ 30 mg/ml	1 : 1	>99%
Human Serum	1 : 1	>90%
Human Serum, 1% Surfactant	1 : 1	>95%
Nucleic Acid Recovery	ProCipitate™: Sample	Recovery
Calf Thymus DNA, $A_{260} = 1.00$	1 : 1	>95%
Total RNA, $A_{260} = 1.00$	1 : 1	>99%

Application	Sample Size	ProCipitate™
Whole Blood, GuSCN Direct Lysis	50 $\mu$ l	250 $\mu$ l
Whole Blood, from dried blood card	7 mm punch	250 $\mu$ l
Plasmid DNA, 96 well format	0.3 ml 2xYT culture	20 $\mu$ l
Plasmid or BAC DNA, bulk processing	200 ml LB culture	5 ml
BAC/PAC DNA, miniprep	1.5 ml 2xYT culture	80 $\mu$ l
Mouse Tail Genomic DNA	Mouse Tail	250-500 $\mu$ l
Plant Genomic DNA	4 mm leaf	100-200 $\mu$ l

## PROTOCOL

1. Resuspend ProCipitate™ by shaking well prior to use.
2. Lyse sample to dissociate nucleic acids from histones and other proteins. Using wide bore or cut pipette tips, add the appropriate volume of ProCipitate™ to deproteinize sample. Use Table above as a guide for volume addition or try several volume ratios starting with a maximum of 1 ml ProCipitate™ to 1 ml of the sample (1 : 1 volume ratio).
3. Gently mix by inversion for 5 minutes at room temperature.
4. Centrifuge sample at 3000 x g for 15 minutes or microfuge at 16,000 x g for 5 minutes.
5. Recover purified nucleic acids contained in the supernatant.
6. Continue with alcohol precipitation or other suitable methods. Note: Buffer condition may be at a moderately acidic pH and there may be a small volume dilution.

### Genomic DNA & Total RNA

ProCipitate™ can be generally substituted for phenol/chloroform in protocols recommended for genomic DNA and total RNA isolations. ProCipitate™ maintains its reactivity in solutions containing surfactant (i.e. SDS, Triton, Tween, N-Lauroylsarcosine). It is suitable for removal of Proteinase K and ribonuclease activity. ProCipitate™ is also available as an integral component of ProPrep™ Genomic & Genomic 96 specifically for whole blood.

### Viral Nucleic Acid Isolation

Reference: Schwab, K.J., De Leon, R., and Sobsey, M.D., *Concentration And Purification Of Beef Extract Mock Eluates From Water Samples For The Detection Of Enteroviruses*,



BIOTECH SUPPORT GROUP

*Hepatitis A Virus, and Norwalk Virus by Reverse Transcription-PCR*, Applied and Env. Microbio, 61:531-537, 1995.

Another polyelectrolyte reagent, Viraffinity™, also can be utilized. Please contact us.

#### ProCipitate™ is Scalable

The volumetric ratio of ProCipitate™ to sample can be adjusted up or down depending on the concentration of protein in the sample. Please contact us to obtain guidance in achieving the optimal volume ratio.

#### Mouse Tail DNA

This protocol comes from email communication with Dr. Tom Nugent at the Scripps Research Institute in San Diego, CA: "We have routinely used this protocol in our laboratory for the last 3 years because it's a relatively simple and very reliable way to process a large number of samples."

1. Digest mouse tail tip in 0.6 ml SET buffer [1% w/v SDS, 5 mM EDTA, 10mM TRIS, 100 mM NaCl] with Proteinase K to a final concentration of 200 µg/ml (2-3 h at 50°C or overnight at 37°C). Spin down tubes. Transfer 0.5 ml into a new tube.
2. Add 0.5 ml ProCipitate™. Incubate 10 min. Mixing every 2-3 min.
3. Spin down sample. Remove 0.5 ml of supernatant.
4. Precipitate DNA with 1/10 volume NaOAc and 0.8 ml Ethanol. Spin down.
5. Dry residual ethanol and resuspend in water or TE.

## References

1. U.S. Patent Numbers 5,294,681, 5,453,493 and other patents pending.
2. U.S. Patent Number 5,538,870, [Method for Preparing Nucleic Acids For Analysis And Kits Useful Therefore](#).
3. [Transgenic labeling of hair cells in the zebrafish acousticolateralis system Brian M. McDermott Jr Gene Expression Patterns Volume 10, Issues 2-3](#), February-March 2010, Pages 113-118
4. Formation of Deoxyguanosine Cross-Links from Calf Thymus DNA Treated with Acrolein and 4-Hydroxy-2-nonenal. Ivan D. Kozekov, Robert J. Turesky, Guillermo R. Alas, Constance M. Harris, Thomas M. Harris, Carmelo J. Rizzo Chemical Research in Toxicology 2010 23 (11), 1701-1713
5. [Genome-wide sequence and functional analysis of early replicating DNA in normal human fibroblasts](#) Stephanie M Cohen, Terrence S Furey, Norman A Doggett, and David G Kaufman BMC Genomics. 2006; 7: 301. Published online 2006 November 29. doi: 10.1186/1471-2164-7-301.
6. Dr. Domon, National Agricultural Research Center for Kyushu Okinawa Region, Japan, Extraction of Rush DNA, unpublished personal correspondence, 2004.
7. Krupey, J., et al, 100,000+ PCRs Possible from 10 ml Blood, poster Biotechniques Symposium, 2003.



BIOTECH SUPPORT GROUP

## CONTACT US

We welcome your questions and comments regarding our products.

Call 732-274-2866, 800-935-0628 (North America) Mon – Fri 9am-6pm EST.

Fax 732-274-2899

Email [sales@biotechsupportgroup.com](mailto:sales@biotechsupportgroup.com)

Mail 1 Deer Park Drive, Suite M, Monmouth Junction, NJ 08852