Minute™ Yeast Mitochondria Enrichment Kit

Catalog number: YM-017

Description

Traditional protocols for yeast mitochondria isolation/extraction involve a range of centrifugation-based subcellular fractionation procedures. Typically the techniques include spheroplast preparation, glass-bead lysis using a homogenization instrument, differential centrifugation and several density gradient procedures using a variety of gradient media with ultracentrifugation. The procedures are very tedious and time consuming. Here we feature a simple and rapid protocol for yeast mitochondria enrichment. The procedure is gentle and instrument-free. Nativ mitochondrial proteins can be isolated from yeast in about one hour without ultracentrifugation. This kit contains optimized detergent-free protein extraction buffers. The protein yield is in the range of 150-250 μ g/sample. The materials provided are sufficient for 50 extractions.

Applications

Proteins extracted with this kit can be used for many downstream applications such as SDS-PAGE analysis, Western blotting, IP. ELISA, enzyme activity assays, proteomics and other biochemical analysis.

Kit components

- 1. 30 ml buffer A
- 2. 10 ml buffer B
- 3. 5 g protein extraction powder
- 4. 4 pestles for 1.5 ml microcentrifuge tube
- 5. 1.5 ml microcentrifuge tube X 50

Storage: Store the kit at -20°C.

Additional Materials Required

Table-Top Microcentrifuge with a maximum rpm of 14,000-16,000. 1 X PBS.

Important Product Information

Prior to plasma membrane isolation addition of protease inhibitor cocktail to buffer A is recommended. For determination of protein concentration, BCA kit (Pierce) is recommended. To study protein phosphorylation, **phosphatase inhibitors** (such as PhosStop from Roche) should be added to buffer A prior to use.

Protocol

- 1. Harvest yeast cells (*S. cerevisiae or S. pombe*) in log growth phase by centrifugation. Collect yeast cells in a 1.5 ml microfuge tube provided. Make sure that the wet volume of pellet is between 30-40 μl. The volume can be easily estimated by comparing it to a 1.5 ml tube with 30-40 μl water.
- 2. Resuspend the pellet in 1 ml cold water and add 100 mg protein extraction powder to the tube. Vortex the tube briefly and centrifuge at top speed in a microcentrfuge for 2 min. **Remove supernatant completely**.
- 3. Grinding the pellet repeatedly with the pestle provided for about 2 min with twisting force. Add 300 µl buffer A to the tube and continue to grind for about thirty seconds (note: The pestle is reusable, for cleaning simply soak it in bleach, rinse with water and dry it with paper towel). Cap the tube and vortex vigorously for 10 seconds.
- 4. Centrifuge the tube at 5000 rpm for 2 min at 4°C. Transfer the supernatant to a fresh pre-chilled microfuge tube and place on ice. Repeat step 3 one more time and combine the supernatants (600 μl) in the microfuge tube on ice.
- 5. Centrifuge the tube at top speed for 20 min at 4°C. Transfer the supernatant to a fresh microfuge tube (this is yeast cytosolic protein fraction). Resuspend the pellet in 200 μl buffer B by vortexing vigorously. Centrifuge the tube at 8000 rpm for 5 min at 4°C.
- 6. Transfer the supernatant to a pre-chilled 1.5 ml tube used in your Lab and add 1 ml cold 1 X PBS to the tube. Invert the tube a few times and centrifuge at top speed at 4°C for 30 min. The pellet contains enriched yeast mitochondria. It can be dissolved in a detergent of your choice depending upon downstream applications.

Application tips: The final protein yield is proportional to grinding frequency and time in step 3. The pestles fit the best with 1.5 ml tubes provided.