

### Oncology and Cytogenetic Products

OaCP S.R.L.



Viale Fanin, 48, Bologna

m@il: ask@oacp.it

www.oacp.it

The DoMo Genetics Line



Improve the efficiency of ISH test:

reduce the costs and save time



## **REF:** 031114-DO

# **Rapid-ISH Integra**

**Contains:** 

1 Vial 100ul (20 test)

All the reagents included in the DoMo Genetics line are protected by current intellectual property laws.

Patent pending (Di Oto E., Monti V., Asioli S.)

#### FOR THE MAXIMUM PERFORMANCE CHANGE YOUR REAGENT BEFORE THE PROTOCOL START

Place the slides in a dry oven at 65 ° C for 30 minutes

- Place in a dry oven at 65  $^{\rm o}$  C a coplin jar with 50ml of Xylene
- Pre reheated a coplin jar with 100ml of SSC2X in the water bath at 77  $^{\circ}$  C
- Pre reheated a coplin jar with 100ml of SSC2X in the water bath at 47  $^{\circ}$  C
- Pre heat in the water bath a coplin jar with 100ml of SSC2X / 1.5% NP40 at 75  $^{\circ}$  C
- On the hybridization plate set the fixed temperature at 75  $^{\circ}\,\mathrm{C}$

Incubate the slides on the plate at 75 ° C for 5 minutes,

- Immerse the slides in Xylene in oven at 65  $^{\rm o}$  C for 30 minutes
- Proceed with 3 sequential washings of the slides with 50ml of xylene in a coplin at RT for 3 minutes / cad.

Dry the slides at RT for 5 minutes

Dehydrate the slides in 2 sequential steps in coplin with 50 ml of 100% ethanol for 5 minutes / Cad.

Dry the slides at RT for 5 minutes

Incubate the slides in Coplin with 2X SSC at 77 °C for a time to be determined between 12 and 18 minutes in relation to the characteristics of the sample

Dissolve 630ul of Proteinase K in the coplin with SSC at 47  $^{\circ}\,\text{C}$ 

- Incubate the slides in the coplin at 47 ° C for a time to be determined between 12 and 18 minutes in relation to the characteristics of the sample
- Then wash the slides in a quick dip into a coplin with 50 ml of  $\ensuremath{\mathsf{SSC2X}}$
- Dehydrate the slides in 3 sequential steps in a coplin with 50 ml of Ethanol 70\% -85\% -100\% for 1 minute / Cad.

Dry the slides at RT for 5 minutes

### On each slide affix 3 ul of probe and 5ul of **Rapid-ISH Integra Buffer**

Cover the area with a cover slip and seal with rubber cement

Set on the hybridization plate a protocol which provides: Denaturation, temperature and time according to the prob specifications; Hybridization, temperature according to the probe specifications, *time : 40 minutes* 

Remove the coverslip and quickly wash slides in a Coplin with 50 ml of SSC2X at RT

Dip the slides in the coplin with SSC2X / 1.5% NP40 at 75  $^{\circ}$  C for 3 minutes

quickly wash slides in a coplin with 50 ml of SSC2X at RT

Dehydrate the slides in 3 sequential steps in a coplin with 50 ml of Ethanol 70% -85% -100% for 1 minute / Cad.

Dry the slides at RT for 5 minutes Affix 5-10 ul of DAPI on each slide, cover with coverslip

Ready for the observation under the microscope

Note

The **RapidISH Integra** buffers should be brought to a temperature of 37  $^{\circ}$  C before use and resuspended to allow optimal mixing of the components.

Solutions should be stored between + 4  $^\circ$  C and -20  $^\circ$  C for better durability.

The transport of the solutions can take place at RT and / or + 4  $^{\circ}$  C.



BUFFER Rapid-ISH Integra, Her2 gene amplification