

Ⅳ

SeleCTEV[™] DNA Enrichment Kit

User Manual – Version 2020/12 - Rev.02



24 reactions



EXO-SEL-LV

EXOSOMICS SpA

Address: Strada del Petriccio e Belriguardo 35, 53100 Siena, Italy

Tel: +39 0577 381308/ Fax +39 0577 381452

Business Registry: SI-137798 VAT#: IT01317090528

info@exosomics.eu www.exosomics.eu

INDEX

Intended Use	3
Product Description	3
Materials Provided/Required	3
Kit components	3
Materials Required but Not Provided	4
Method Description and procedure	4
Method Description	4
Procedure	5
Troubleshooting	7
Warnings and Precautions	8
Safety information	8
General precautions	8
Limitations	8
Symbols	8
Contact Information	9

Intended Use

SeleCTEV[™] DNA Enrichment Kit Enrichment kit is a pre-analytical kit to purify both circulating cell free DNA and tumor-originated DNA from tumor enriched extracellular vesicles (EVs) and exosomes from plasma.

The **SeleCTEV[™] DNA Enrichment Kit** (EXO-SEL-LV) kit is optimized for input volumes ranging from 0.5 ml up to 2 ml of plasma.

The kit is designed to be used with any downstream application employing enzymatic amplification or enzymatic modifications of DNA isolated with SeleCTEV, followed by signal detection or amplification (such as PCR, real time PCR and digital PCR).

Exosomics **SeleCTEV**TM **DNA Enrichment Kit** does not provide a diagnostic result. It is responsibility of the user to use and validate the kit in conjunction with any downstream assay.

Product Description

SeleCTEVTM **DNA Enrichment Kit** is ready-to-use and it is meant for running 24 tests. Kit allows the selective isolation of tumor-originated nucleic acids from tumor enriched EVs and exosomes, from a minimum of 500 μ l of plasma following two subsequent working steps:

- 1. Cell free circulating and EV-associated DNA isolation from biofluids.
- 2. DNA purification.

The purification is based on Exosomics proprietary peptide affinity method and does not require any special equipment, such as ultracentrifugation or chromatography, with a turnaround time of a minimum of 4 hours.

Materials Provided/Required

Kit components, meant to run a total of 24 reactions, and their storage conditions are listed below table.

Kit components

Product Code (Input Volume)			EXO-SEL-LV (0,5-2ml)
Components	Name	Description	Quantity (volume/number)
Isolation Agent	EXO-IA	Lyophilized Reagent for isolation	1 vial (2.5 mg)
Resuspension Buffer	EXO-RB	EXO-IA Resuspension Buffer	1 vial (1 ml)
10X Isolation Buffer	EXO-IB	Diluent for Isolation	1 bottle (30 ml)
Isolation Tubes	EXO-IsoT-2ml	Tubes for EV isolation	48 tubes (2 ml)
Proteinase K	EXO-PK	Reagent for Protein Digestion	1 ml (20 mg/ml)
Lysis Buffer	EXO-LB	Buffer for vesicle Lysis	1 bottle (10 ml)
Washing Buffer 1	EXO-WB1	Buffer for column Washing	1 bottle (15 ml)
Washing Buffer 2	EXO-WB2	Buffer for column Washing	1 bottle (15 ml)
Elution Buffer	EXO-EB	Buffer for DNA elution	1 bottle (4 ml)
DNA purification columns	EXO-DC	Columns for DNA Elution	24 columns
Elution Tubes	EXO-ColT-1.5 ml	Tubes for pure DNA collection	24 collection tubes (1.5 ml)

SeleCTEV™ DNA Enrichment Kit Handbook Cat N: EXO-SEL-LV Store all the supplied reagents according to the instructions on the respective packages. In particular:

- Upon EXO-IA resuspension in EXO-RB, store the vial at +2/+8 °C;
- Upon preparation of EXO-WB1 and EXO-WB2 reagents, store these buffers at +2/+8 °C;
- Upon arrival, store EXO-PK at +2/+8 °C;
- If properly stored, all the reagents provided with the kit are stable until the expiration date printed on the label.

Materials Required but Not Provided

Ш	Protease inhibitor (Example: Sigma cat num. P8340)
	96-100% Ethanol (Example: Sigma Aldrich, Cat. Num.02860) for dilution of EXO-WB1 and EXO-WB2
	buffers and for DNA purification
	Disposable Gloves
	Single-use and/or pipettes with disposable tips
	Pipettes for reagent preparation
	Ultrapure water for dilution of 10X EXO-IB
	Heating block, or water bath for incubation at 56 °C
	Benchtop centrifuge with rotor for 2 ml reaction tubes (kit validated on Eppendorf 5415R) for Low
	Volume kit (EXO-SEL-LV).
	Vortex

Method Description and procedure

Method Description

Method: SeleCTEV[™] DNA kit method is based on Exosomics' proprietary peptide-affinity method that selectively binds and enriches for tumor derived exosomes from which DNA is then extracted. SeleCTEV[™] DNA affinity method also recovers circulating nucleic acids, making it the method of choice to recover the total circulating nucleic acid plus the tumor derived fraction from exosomes.

SeleCTEV[™] DNA Enrichment Kit has been optimized for sample volumes ranging from 0.5 ml to 2 ml of plasma. Follow steps 1-4 up to 1 ml of plasma. For volume of plasma >1 ml, the best performance is obtained by splitting plasma 1:1 into two isolation tubes (EXO-IsoT-2 ml) and then proceeding through steps 1-4 as for 1 ml samples.

Sample collection: Exosomics S.p.A. suggests the following recommendations,

- tubes for blood collection with anticoagulant, K2 EDTA, or CTAD or citrate;
- store the collection tubes at room temperature (20-25°C);
- process the samples within 30 minutes or in a shorter time possible after collection;
- centrifuge the blood samples at 1500 g 10 minutes at room temperature (20-25°C);
- recover the plasma on the top, being careful to not contaminate the samples with red blood cells;
- aliquoting is recommended since freeze-and-thaw cycles reduce the quality of the sample;
- store the plasma collected freeze (preferred temperature -80°C) or proceed with the **SeleCTEV**[™] **DNA Enrichment Kit** protocol.

Procedure

SeleCTEV[™] **DNA** allows the isolation of extracellular vesicles and tumor originated exosomes from plasma, and extraction of their DNA through a four-step procedure:

- 1. Reagent preparation
- 2. Plasma preparation
- 3. Cell free DNA and EV isolation from plasma
- 4. DNA purification

1 Reagent preparation:

- 1.2 **Isolation Agent (EXO-IA):** add 1 ml of Resuspension Buffer (EXO-RB) into the Isolation agent (EXO-IA) vial. Gently tap the vial and visually check for complete resuspension of the lyophilized reagent. Do not pipet up and down.
- 1.3 **1X Isolation Buffer (1X-IB):** dilute **10X Isolation Buffer (EXO-IB)** in fresh ultrapure water to a final 1X concentration (i.e. 1 ml of EXO-IB and 9 ml of ultrapure water) and label the vial as "1X-IB".
- 1.4 Washing Buffer 1 (EXO-WB1): add 9.4 ml of pure Ethanol (96-100%) to EXO-WB1 bottle (15 ml). Mix well by inverting 6-8 times.
- 1.5 **Washing Buffer 2 (EXO-WB2):** add 10.5 ml of pure Ethanol (96-100%) to EXO-WB2 bottle (15 ml). Mix well by inverting 6-8 times.

2 Plasma preparation:

- 2.1 Pre-clear plasma sample by centrifuging at 1200 g for 20 min at 10 °C to eliminate red blood cells and cellular debris.
- 2.2 Discard the pellet and debris and transfer the supernatant in the appropriate tube (EXO-IsoT-2 ml, 2 ml volume tubes for isolation with SeleCTEV[™] DNA Enrichment Kit)
- 2.3 Dilute pre-cleared plasma with **1X-IB** buffer to a final 1:1 dilution (i.e. If processing 0.5 ml of plasma, add 0.5 ml of 1X-IB). If sample volume is >1 ml, split the sample into two isolation tubes (EXO-IsoT-2 ml) and then proceed to a final 1:1 dilution with **1X-IB** (i.e. If processing 2 ml of plasma, split sample equally in two isolation tubes and then add 1 ml of 1X- IB to each tubes).
- 2.4 Add protease inhibitor cocktail to each sample (1:1000 v/v protease: diluted plasma. Not provided with the kit, we recommend Sigma cat num. P8340.)

3 Cell-free DNA and EV isolation from plasma:

3.1. Add resuspended Isolation agent (EXO-IA) to each vial of pre-cleared diluted sample.

EXO-SEL-LV: add 20 μl of EXO-IA to pre-cleared diluted plasma

- 3.2. Mix well by pipetting and inverting the tube.
- 3.3. Incubation time is 2 hours at RT under rotation.
- 3.4. Centrifuge 15 min at 16 000 g at RT.
- 3.5. Discard the supernatant, carefully avoiding to dislodge the pellet. Eliminate the remaining supernatant from the tube with a pipette.
- 3.6. Gently add 1 ml of **1X Isolation Buffer (1X-IB)** directly on the pellet, without disrupting it. Spin the sample at 7000 g for 7 min at RT.
- (\mathbf{i}) If the pellet is not visible at this step, refer to Technical Support (Refer to Troubleshooting, p.6).
- 3.7. Repeat steps 3.5-3.6 one more time.

- 3.8. Discard the supernatant, carefully avoiding to dislodge the pellet. Eliminate the remaining supernatant from the tube with a pipette.
- 3.9. Resuspend (each) pellet in 200 µl of Isolation Buffer (1X-IB).
- (i) We advise to proceed directly to step 4 (DNA purification) to obtain optimal DNA recovery.

4 DNA purification:

4.1 EV Lysis:

- 4.1.1 Add 20 μl of Proteinase K (20 mg/ml), **(EXO-PK)**, to each resuspended pellet and mix by gently vortexing the tube.
- 4.1.2 Add 200 µl of Lysis Buffer (EXO-LB) to each tube
- (i) If processing >1 ml of plasma, add 200 μl of EXO-LB to each isolation tube (EXO-IsoT-2 ml)
- 4.1.3 Mix well by vortexing 30 sec.
- 4.1.4 Incubate samples at 56 °C for 1 hour.

4.2 DNA purification:

- 4.2.1 Add 200 μl of Ethanol 96-100% to each tube and mix by briefly vortexing the tube.
- 4.2.2 Transfer the mixtures in a **DNA Spin Column (EXO-DC)** and centrifuge at 10 000 g for 1 min.
- If processing >1 ml of plasma, repeat steps from 4.1.1 to 4.2.1 with two tubes, and then load them into the same DNA Spin Column (4.2.2).
- 4.2.3 Discard the flow-through.
- 4.2.4 Add 500 μ l of **Washing Buffer 1 (EXO-WB1)**, centrifuge at 10 000 g for 1 min and discard the flow-through.
- 4.2.5 Add 500 μ l of **Washing Buffer 2 (EXO-WB2)**, centrifuge at 10 000 g for 1 min and discard the flow-through.
- 4.2.6 Centrifuge 2 additional min at 16 000 g.
- 4.2.7 Transfer the spin column to an **Elution Tube (EXO-ColT-1.5 ml)**.
- 4.2.8 Elute the DNA from the column adding 20 μl of **Elution Buffer (EXO-EB)**.
- 4.2.9 Incubate for 5 min at RT.
- 4.2.10 Centrifuge 1 min at 14 000 g.
- 4.2.11 Transfer the eluate to the spin column, incubate for 5 minutes.
- 4.2.12 Repeat point 4.2.10 one more time.
- 4.2.13 Samples can now be used for further analyses or stored at -20 °C.

SeleCTEV™ DNA Enrichment Kit Handbook Cat N: EXO-SEL-LV

Troubleshooting

This table may solve some technical problems that could arise during **SeleCTEV**TM **DNA Enrichment Kit** protocol execution. For more information, please contact us at support@exosomics.eu.

Technical	Potential Causes	Suggestions and comments
Problems		
	Poor plasma quality due to delayed blood processing	Repeat blood processing according to Sample processing in "Method Description and Procedure" section.
	Plasma samples are frozen and thawed multiple times	Always use fresh samples or samples thawed once.
	Prolonged sample storage at room temperature	Do not keep the samples at RT for prolonged time.
	Incomplete resuspension of the peptide	Peptide solution may initially look cloudy after resuspension in resuspension buffer (EXO-RB). Do not vortex the solution, simply tap the vial to resuspend the peptide. Make sure that the peptide is fully resuspended in EXO-RB and the final solution looks clear.
Poor DNA	No visible pellet	It may occasionally occur but should not affect DNA recovery.
recovery	Lysis buffer (EXO-LB) and pellet- proteinase K solutions not sufficiently mixed	Mix lysis buffer (EXO-LB) and pellet-proteinase K solution well by pipetting up and down and vortexing at least 30" to completely resuspend the peptide pellet.
	Inefficient sample lysis	Use fresh proteinase K. If needed, increase incubation time with proteinase K.
	Sub-optimal ethanol percentage	Use fresh 96-100% ethanol. Do not use denatured alcohol which may contain methanol.
	Clogged DNA spin column	Repeat the procedure increasing the incubation time in proteinase K.
	Wash buffers 1 and 2 (EXO-WB1; EXO-WB2) prepared incorrectly	Check that these buffers were diluted in the correct volume of 96-100% ethanol (see page 7).
	The eluate volume is lower than the applied volume	Expect to recover an eluate volume with 2-3 μ l less than the applied volume due to retention of the silica membrane.
	Presence of ethanol traces in eluate	Make sure to remove all ethanol residuals from the column (EXO-DC) before eluting the sample.
	Extremely low or no DNA recovered	See poor DNA recovery section above for troubleshooting.
DNA not suitable for	Not optimized elution volume	Calculate the optimal elution volume for PCR reaction.
enzymatic reaction	New PCR assay	If the PCR assay is changed, readjust the eluate volume.
	Interference due to plasma inhibitors	Consider the presence of plasma inhibitors such as natural or synthetic small molecule (therapeutics) that may end up in the eluate and inhibit DNA amplification.

Warnings and Precautions

SeleCTEV[™] **DNA Enrichment kit** does not provide a diagnostic result. It is the sole responsibility of the user to use and validate the kit in conjunction with a downstream assay.

All products sold by Exosomics are subject to extensive quality control procedures and are warranted to perform as described and used correctly. Any problems should be reported immediately. The kit contents are for laboratory use only, and they must be stored in the laboratory and must not be used for purposes other than intended.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles.

Lysis Buffer (EXO-LB) and Washing Buffer 1 (EXO-WB1) contain guanidinium salts, and should be handled with care. Guanidinium salts form highly reactive compounds when combined with bleach, thus care must be taken to properly dispose of any of these solutions.

For more details, please refer to the product **SeleCTEV[™] DNA Enrichment kit** Safety Data Sheet.

General precautions

Always operate in accordance with Good Laboratory Practice (GLP) guidelines and the instructions included in this Handbook.

Handle and dispose of waste materials and reagents from the use of this kit in accordance with current regulations on staff safety and environmental protection.

Limitations

Dilution of the reagents, other than as described in this handbook, is not recommended and will result in a loss of performance.

Attention should be paid to expiration dates and storage conditions printed on the box and labels of all components. Do not use expired or incorrectly stored components.

Symbols

Σ	Contains reagents sufficient for n tests
REF	Catalogue number
LOT	Batch or lot code
QUAN	Quantity
VOL	Volume
	Use by
X	Temperature limit
[]i	Operating instructions
	Manufacturer
C€	CE Mark

Contact Information

Exosomics S.p.A.

Address: Strada del Petriccio e Belriguardo 35, 53100 Siena, Italy

Business Registry: SI-137798 **VAT#:** IT01317090528

Administration: Tel: +39 0577 381308

For info: administration@exosomics.eu

Technical Support:

For info: support@exosomics.eu

Sales & Marketing: Tel: +39 0577 381309

SeleCTEV™ DNA Enrichment Kit Handbook Cat N: EXO-SEL-LV