



11 Park Drive, Suite 12
Boston, MA 02215

Hamster monoclonal anti mouse CD152 (CTLA-4) (Azide-free & Low endotoxin)

ORDERING INFORMATION

Catalog Number:	gAP-0028
Size:	1.00 mg
Storage:	< -20° C
Immunogen:	Mouse CTLA-4 protein
Ig Type:	Hamster IgG
Clone	AP-MAB0835
Endotoxin Level	< 0.002EU/μg IgG*
Applications:	FC and Blocking

Description: CD152 is a 33 kD member of the immunoglobulin superfamily also known as CTLA-4 or Ly-56. It is expressed on activated T and B lymphocytes. CD152 is similar to CD28 in amino acid sequence, structure, and genomic organization and these two receptors share common B7 family counter-receptors (B7-1, B7-2). Whereas CD28 delivers a costimulatory signal in T cell activation, CTLA-4 negatively regulates cell-mediated immune responses. CD152 is thought to play a role in the induction and maintenance of immunological tolerance, as well as the development of protective immunity and thymocyte regulation.

Preparation: This antibody was produced from a hybridoma (mouse myeloma fused with spleen cells from a Hamster immunized with **Mouse CTLA-4 protein**).

Formulation: The IgG fraction of **culture supernatant** was purified by Protein A/G affinity chromatography and lyophilized from a 0.2 μm filtered solution in phosphate-buffered saline (PBS, **Azide Free**).

Reconstitution: Reconstitute the antibody with sterile PBS and the reconstituted antibody can be aliquoted and stored frozen at < -20 for at least for six months without detectable loss of activity. **Avoid repeated freeze-thaw cycles**. Lyophilized samples are stable for 2 years from date of receipt when stored at -70°C.

***Endotoxin Level:** **Extremely low level of LPS (< 0.002EU/μg IgG)**

Application(s):

- 1. FC**
- 2. In vitro stimulation**
- 3. In vitro and in vivo blocking of ligand binding**

*** The antibody is produced by in vitro culture.**

Contact & Ordering Information: Angio-Proteomie, 11 Park Drive, Suite 12, Boston, MA 02215, USA. Tel: 617-549-2665; Fax: (480) 247-4337, angioproteomie@gmail.com