



11 Park Drive, Suite 12
Boston, MA 02215

Rat anti-human CCR3

ORDERING INFORMATION

Catalog No: hAP-0063
Size: 100 µg
Storage: < -20° C
Immunogen: Human CCR3 transfected 293 cells
Ig Type: Rat IgG2
Clone: MM0063-7F12
Applications: FC and IHC (Frozen)

Preparation: This antibody was produced from a hybridoma (mouse myeloma fused with spleen cells from a Rat immunized with 293 cells transfected with human CCR3).

Formulation and Storage: The IgG fraction of culture supernatant was purified by Protein G affinity chromatography and lyophilized from a 0.2 µm filtered solution in phosphate-buffered saline (PBS). Lyophilized samples are stable for 2 years from date of receipt when stored at -70°C.

Reconstitution: Reconstitute the antibody with 200 µl sterile PBS and the final concentration is 500µg/ml. Reconstituted antibody can also be aliquoted and stored frozen at < -20 for at least for six months without detectable loss of activity. **Avoid repeated freeze-thaw cycles.**

Specificity: This antibody was selected for its ability to detect human CCR3, not other CCRs.

Applications:

- | | | | |
|----|---------------------|------------|--------------------|
| 1. | FC: | Yes | (1:100-800) |
| 2. | IHC (Frozen) | Yes | (1:50-200) |

Final dilutions should be optimized by the end users and the dilutions provided in the data sheet are for reference only

Contact & Ordering Information: Angio-Proteomie, 11 Park Drive, Suite 12, Boston, MA 02215, USA. Tel: 001-6175492665; Fax: 001-4802474337, angioproteomie@gmail.com