



11 Park Drive, Suite 12  
Boston, MA 02215

# Endothelial Growth Medium

## ORDERING INFORMATION

**Product Name**            **Endothelial Growth Medium (ENDO-Growth Medium)**  
**Catalog No:**            **cAP-02**  
**Size:**                      **500ml**  
**Storage:**                 **4°C**

### Description

Endothelial growth medium (ENDO-Growth medium) is freshly prepared by mixing the two components: (1) 475ml of ENDO-Basal medium (cAP-03) and (2) 25ml of ENDO-Growth supplement (cAP-04), which contains the essential nutrients for endothelial cell proliferation (including 5% fetal bovine serum, recombinant growth factors, and 1x penicillin and streptomycin).

### Shipping Condition:

Ambient temperature (blue ice, seasonally)

**Storage Condition:** The ENDO-Growth medium is stored at 4°C. A change in color or appearance of precipitate may indicate deterioration.

**Shelf Life:** The ENDO-Growth medium: at least 1 month from the date of receipt under proper storage condition.

### Use of ENDO-Growth Medium for culturing human endothelial cells:

1. Take the amount of ENDO-growth medium needed and warm to 37°C in a water bath or incubator.
2. Avoid frequent temperature change to the entire bottle of medium.
3. The ENDO-growth medium is stable for at least 1 month stored at 4°C. We recommend changing medium every 2 days for regular cell culture.

## THIS PRODUCT IS FOR RESEARCH USE ONLY NOT FOR DIAGNOSTIC OR THERAPEUTIC APPLICATIONS

**Caution:** Handling human derived products is potentially biohazardous. Although each cell strain testes negative for HIV, HBV and HCV DNA, diagnostic tests are not necessarily 100% accurate, therefore, proper precautions must be taken to avoid inadvertent exposure. Always wear gloves and safety glasses when working these materials. Never mouth pipette. We recommend following the universal procedures for handling products of human origin as the minimum precaution against contamination.

Contact & Ordering Information: Angio-Proteomie, 11 Park Drive, Suite 12, Boston, MA 02215, USA. Tel: 001-6175492665; Fax: 001-4802474337, [angioproteomie@gmail.com](mailto:angioproteomie@gmail.com)