

Certificate of Analysis

Product:	Rabbit Monoclonal Antibody Anti-SARS-CoV-2 (2019-nCov) Nucleocapsid Rabbit Monoclonal Antibody, Clone RM420
Catalog No.:	31-1307-00
Lot No.:	
Clone:	RM420
Isotype:	Rabbit IgG
Specificity:	This antibody reacts 2019-nCov Nucleocapsid Protein.
Application:	Sandwich ELISA assay as a capture antibody, pair with chimeric human monoclonal anti-2019-nCov NP clone RMH03
Immunogen:	A recombinant 2019-nCov Nucleocapsid Protein expressed in <i>E.coli</i> .
Purity:	Protein A affinity purified from an animal origin-free culture supernatant.
Size:	
Concentration:	1 mg/mL
Buffer:	PBS Only
Storage and Stability:	Store vial at 4° C for short term. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. This product is stable for 1-2 weeks at 4° C as an undiluted liquid.
Country of Origin:	U.S.A.
Intended Use:	For Research Use Only Not for Diagnostic or Therapeutic Use

RELATED PRODUCT:

Anti-SARS-CoV-2 (2019-nCov) Nucleocapsid Chimeric Human Monoclonal Antibody, Clone RMH03