

Survivin 2B-Specific Antibody

CH3 Catalog # 150001, Unit Size: 100 µl

- PRODUCT OVERVIEW** : Anti-Survivin 2B is an IgY-purified chicken polyclonal antibody that was raised against survivin 2B peptide (CZRITREEHKKHSS).
- PURIFICATION** : Affinity purification by agarose chromatography using conjugated immunizing peptide.
- FORMULATION** : PBS/glycerol/0.02% sodium azide.
- RECOMMENDED USAGE** : Western Blotting.
- TESTING** : Tested using protein extracts from HeLa cells transfected with either survivin, survivin 2B or survivin ΔEx3. This reagent recognizes only survivin 2B. For experimental procedures, consult Li, F. and Ling, X. 2006 (<http://www.ncbi.nlm.nih.gov/pubmed/16557517>).
- STORAGE** : Avoid freeze/thaw cycles. Store at 4°C or prepare appropriately sized aliquots and keep at -20°C.

All products are for Research Use Only and not for in vitro diagnostic use

This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.

University Commons • Suite 8 • 1416 Sweet Home Road • Amherst, NY 14228 • 716.688.5222
www.ch3biosystems.com



삼보메디칼(www.sambomed.co.kr)

Tel: 02) 575-4945, Fax : 02) 574-7595