

Chicken Anti-Survivin and Anti-2B-Specific(but not DEx3) IgY Antibody

CH3 Catalog # 150003, Unit Size: 100 µl

- PRODUCT OVERVIEW : NA
- PURIFICATION : Affinity purification by agarose chromatography using conjugated immunizing peptide.
- FORMULATION : PBS/glycerol/0.02% sodium azide.
- RECOMMENDED USAGE : Western Blotting.
- TESTING : Tested using protein extracts from HeLa cells transfected either with survivin, survivin 2B or survivin ΔEx3. This reagent recognizes only survivin 2B. For experimental procedures, consult Li, F. and Ling, X. 2006 (<http://www.ncbi.nlm.nih.gov/pubmed/16557517>).
- STORAGE : Avoid freeze/thaw cycles. Store at 4°C or prepare appropriately sized aliquots and keep at -20°C.

*All products are for Research Use Only and not for in vitro diagnostic use
This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.*

University Commons • Suite 8 • 1416 Sweet Home Road • Amherst, NY 14228 • 716.688.5222
www.ch3biosystems.com



삼보메디칼(www.sambomed.co.kr)

Tel: 02) 575-4945, Fax : 02) 574-7595