

HepG2 Cell Line

CH3 Catalog # 930001, Unit Size: Flask or Vial

PRODUCT OVERVIEW	: Human hepatocellular carcinoma cell line.
PURIFICATION	: NA
FORMULATION	: 25 to 50 X 10 ⁶ cells in culture flask containing RPMI 1640 culture medium and 5% fetal bovine serum or frozen cell pellet in cryo-protectant.
RECOMMENDED USAGE	: Cell culture. Useful to provide cell proteins for large-scale purification in downstream applications. Cells adapt well to a wide variety of culture media. Fetal bovine serum at 5% is recommended.
TESTING	: Cultured on tissue culture plastic in 5% fetal bovine serum with either DMEM or RPMI 1640 culture media.
STORAGE	: Place in culture immediately upon arrival: 35° humidified CO ₂ incubator. Use 5 to 7.5% CO ₂

For Research Use Only. Not for in vitro diagnostic use

University Commons • Suite 8 • 1416 Sweet Home Road • Amherst, NY 14228 • 716.688.5222
www.ch3biosystems.com



Sambo

삼보메디칼(www.sambomed.co.kr)

Tel: 02) 575-4945, Fax : 02) 574-7595