

PC-3 Cell Line

CH3 Catalog # 940012, Unit Size: Flask or Vial

PRODUCT OVERVIEW	: Human prostate cancer cell line.
PURIFICATION	: NA
FORMULATION	: 10 to 20 X 10 ⁶ cells in flask containing RPMI 1640 culture medium and 5% fetal bovine serum or frozen cell pellet in cryoprotectant.
RECOMMENDED USAGE	: Cell culture. Useful for large-scale cell culture to provide human cellular proteins for purification and analysis in downstream applications. Adapt well to a wide variety of culture media. Fetal bovine serum at 5% is recommended.
TESTING	: Cultured on tissue culture plastic in 5% fetal bovine serum with either DMEM or RPMI 1640.
STORAGE	: Place in culture immediately upon arrival: 35-37° humidified CO ₂ incubator. Use 5 to 7.5% CO ₂ .

*All products are for Research Use Only and not for in vitro diagnostic use
This Package Insert is provided for product evaluation purposes only and is not intended to be
used in place of the Package Insert shipped with the product.*

University Commons • Suite 8 • 1416 Sweet Home Road • Amherst, NY 14228 • 716.688.5222
www.ch3biosystems.com



삼보메디칼(www.sambomed.co.kr)

Tel: 02) 575-4945, Fax : 02) 574-7595