

Cationic DOTAP/DC-CHOL Liposomes Adjuvants

Code No. : F50101V

Size : 2 ml

Image
Coming
Soon

Product Description

Specifications

Product code: F50101V

Lipid composition: DOTAP/DC-Cholesterol (50:50 mol/mol)

Mean particle size: 50-80 nm (by DSL)

Lipid concentration: 40 mM (24.7 mg/mL)

Bulk solution: D5W (dextrose in milli-Q water)

Shelf-life: 6 months guaranteed for sterility for unopened vials stored at 2-8 degree C
(guaranteed starting from the date of delivery, replacement only)

DOTAP: 1,2-dioleoyl-3-trimethylammonium-propane (chloride salt)

DOPE: DC-Cholesterol • HCl, 3b-[N-(N',N'-dimethylaminoethane)-caramoyl]chelesterol hydrochloride



삼보메디칼(www.sambomed.co.kr)

Tel: 02) 575-4945, Fax : 070) 4104-7595