



## Certificate of Analysis

<b>Product:</b>	Rabbit Monoclonal Antibody Anti-Ki-67 Rabbit Monoclonal Antibody, Clone RM227
<b>Catalog No.:</b>	31-1106-00
<b>Lot No.:</b>	
<b>Clone:</b>	RM227
<b>Specificity:</b>	This antibody reacts to Human Antigen Ki-67.
<b>Application:</b>	Immunohistochemistry
<b>Immunogen:</b>	A peptide corresponding to the internal region of human Ki67
<b>Purity:</b>	Protein A affinity purified from an animal origin-free culture supernatant
<b>Size:</b>	31-1106-00-S: 0.1 mL 31-1106-00-T: 0.25 mL 31-1106-00-R: 1.0 mL
<b>Buffer:</b>	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide
<b>Usage:</b>	IHC: 1:100 dilution.
<b>Storage and Stability:</b>	Stable for 1 Year at -20.0°C from date of receipt
<b>Country of Origin:</b>	U.S.A.
<b>Intended Use:</b>	<b>For Research Use Only Not for Diagnostic or Therapeutic Use</b>



Immunohistochemical staining of formalin fixed and paraffin embedded human breast cancer tissue sections.



Immunohistochemical staining of formalin fixed and paraffin embedded human tonsil tissue sections.