

## Plain 24:0PC/CHOL Liposomes (55:45 mol/mol, 120 nm)

Code No. : F10107

Size : 10 ml

Image  
Coming  
Soon

### Product Description

Plain 24:0PC/CHOL Liposomes (55:45 mol/mol, 120 nm)

This preformed liposome reagent is prepared in 10% sucrose and contains NO drug loading battery. It is typically used as a placebo control for the drug loaded liposomes with the same lipid composition

### Specifications

Product code: F10107

Lipid composition: 24:0PC/CHOL (55:45 mol/mol)

Mean particle size: 110 nm (100-130 nm)

Lipid concentration: 60 mM (60 - 65 mM)

Internal/External solution: 10% sucrose, 20 mM HEPES, pH 7.0 - 7.5, <10% ethanol

24:0PC: 1,2-dilignoceroyl-sn-glycero-3-phosphocholine

CHOL: cholesterol