

Plain Pure DOPC Liposomes (100%, 100 nm)

Code No. : F10109

Size : 10 ml

Image
Coming
Soon

Product Description

Plain Pure DOPC Liposomes (100%, 100 nm)

This preformed liposome reagent is prepared in 10% sucrose and contains NO drug loading battery. It is typically used as a placebo control for the drug loaded liposomes with the same lipid composition

Specifications

Product code: F10109

Lipid composition: DOPC

Mean particle size: 100 nm (90-120 nm)

Lipid concentration: 60 mM (60 - 65 mM)

Internal/External solution: 10% sucrose, 20 mM HEPES, pH 7.0 - 7.5, <10% ethanol

DOPC: 18:1 (Δ^9 -Cis) PC, 1,2-dioleoyl-sn-glycero-3-phosphocholine