

Plain Pure POPC Liposomes (100%, 100 nm)

Code No. : F10110

Size : 10 ml

Image
Coming
Soon

Product Description

Plain Pure POPC Liposomes (100%, 100 nm)

This preformed liposome reagent is prepared in 10% sucrose and contains NO drug loading battery. It is typically used as a placebo control for the drug loaded liposomes with the same lipid composition

Specifications

Product code: F10110

Lipid composition: POPC

Mean particle size: 100 nm (90-120 nm)

Lipid concentration: 60 mM (60 - 65 mM)

Internal/External solution: 10% sucrose, 20 mM HEPES, pH 7.0 - 7.5, <10% ethanol

POPC: 16:0-18:1 PC,1-palmitoyl-2-oleoyl-sn-glycero-3-phosphocholine