

Plain Control Liposomes for Doxoves[®]

Code No. : F30204B-C

Size : 2 ml & 10 ml

Image
Coming
Soon

Product Description

Plain Control Liposomes for Doxoves[®]

This product of plain liposomes contains all the ingredients and has similar physical properties as that for Doxoves[®], except the drug. It contains the same lipid composition and the drug loading battery, ammonium sulfate inside the liposomes.

Specifications :

Product code: F30204B-C

Product name: Control liposomes for Doxoves[®]

Lipid composition: HSPC/CHOL/mPEG2000-DSPE (56.2:38.5:5.3 mol/mol)

Lipid concentration: 40 mM

Dox HCl concentration: none

Loading battery: 250 mM ammonium sulfate

Mean particle size: 85 nm (80-90 nm)

External solution: 10% sucrose, 10 mM histidine, pH 6.0 - 6.5

Storage and Stability: stable for 6 months when stored in refrigerator (2-8 degree C). Use within 1 month weeks after first use or open in BSC to keep sterile for longer stability.

AVOID FREEZING

HSPC: hydrogenated phosphatidylcholine (Soy)

CHOL: cholesterol

mPEG2000-DSPE: 18:0 PEG2000 PE:1,2-distearoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-2000]

Product For Research Only And NOT For Human And Companion Animal Use

제조원 : FormuMax Scientific Inc. (U.S.A)



삼보메디칼(www.sambomed.co.kr)

Tel: 02) 575-4945, Fax : 02) 574-7595