

## Cationic Liposomes DMTAP/DOPE

Code No. : F50101

Size : 2 ml

Image  
Coming  
Soon

### Product Description

#### Specifications

Product code: F50101

Lipid composition: DMTAP/DOPE (50:50 mol/mol)

Mean particle size: 90-110 nm (by DSL)

Lipid concentration: 5 mM (3.5 mg/mL)

Bulk solution: nuclease-free water, 10 vol% ethanol, pH 6.0 - 7.0

Shelf-life: 3 months guaranteed for sterility for unopened vials stored at 2-8 degree C  
(guaranteed starting from the date of delivery, replacement only)

DMTAP: 1,2-dimyristoyl-3-trimethylammonium-propane (chloride salt)

DOPE: 1,2-dioleoyl-sn-glycero-3-phosphoethanolamine



삼보메디칼([www.sambomed.co.kr](http://www.sambomed.co.kr))

Tel: 02) 575-4945, Fax : 070) 4104-7595