

Cationic Liposomes DOTAP/DOPE

Code No. : F50104

Size : 2 ml



Product Description

Specifications

Product code: F50102

Lipid composition: DOTAP/DOPE (50:50 mol/mol)

Mean particle size: 90-110 nm

Lipid concentration: 5 mM (3.61 mg/mL)

Bulk solution: nuclease free water, 10 % ethanol, pH 6.0 - 7.0

Shelf-life: 3 months guaranteed for sterility for unopened vials stored at 2-8 degree C
(guaranteed starting from the date of delivery, replacement only)

DOTAP: 1,2-dioleoyl-3-trimethylammonium-propane (chloride salt)

DOPE: 1,2-dioleoyl-sn-glycero-3-phosphoethanolamine



삼보메디칼(www.sambomed.co.kr)

Tel: 02) 575-4945, Fax : 070) 4104-7595