

Cationic Liposomes DODAP/DOPE

Code No. : F50104

Size : 2 ml



Product Description

Specifications

Product code: F50104

Lipid composition: DODAP/DOPE (50:50 mol/mol)

Mean particle size: 90-110 nm

Lipid concentration: 5 mM (3.48 mg/mL)

Bulk solution: nuclease free water, 10 vol% ethanol, pH 3.0 - 4.0

Shelf-life: 1 month guaranteed for sterility for unopened vials stored at 2-8 degree C

(guaranteed starting from the date of delivery, replacement only)

* Note: DODAP liposomes should be stored at 2-8 degree C at all times and used within one month.

DODAP: 1,2-dioleoyl-3-dimethylammonium-propane

DOPE: 1,2-dioleoyl-sn-glycero-3-phosphoethanolamine



삼보메디칼(www.sambomed.co.kr)

Tel: 02) 575-4945, Fax : 070) 4104-7595