

Cationic Liposomes DC-CHOL/DOPE

Code No. : F50105

Size : 2 ml



Product Description

Specifications

Product code: F50105

Lipid composition: DC-CHOL/DOPE (30:70 mol/mol)

Mean particle size: 90-110 nm (by DSL)

Lipid concentration: 5 mM (3.41 mg/mL)

Bulk solution: nuclease-free water, 10% Ethanol, pH 6.0 - 7.0

Shelf-life: 3 months guaranteed for sterility for unopened vials stored at 2-8 degree C
(guaranteed starting from the date of delivery, replacement only)

DC-CHOL: DC-Cholesterol • HCl, 3b-[N-(N',N'-dimethylaminoethane)-caramoyl]cholesterol hydrochloride

DOPE: 1,2-dioleoyl-sn-glycero-3-phosphoethanolamine



삼보메디칼(www.sambomed.co.kr)

Tel: 02) 575-4945, Fax : 070) 4104-7595