

## Non-Fluorescent DOPC/CHOL Liposomes (300nm)

Code No. : F60103F3-P

Size : 1 ml & 5 ml

### Product Description

#### Placebo Control DOPC/CHOL Liposomes (300nm)

This placebo control liposomes is made to match the lipid composition and characteristics for the non-PEGylated fluorescent liposomes. The product is sterile filtered and filled in sterile vials guaranteed to be stable for 6 months when stored in refrigerator.

#### Specifications

**Lipid composition:** DOPC/CHOL/(54:45 mol/mol)

**Lipid concentration:** 50 mM (50-55 mM)

**Dye :** no dye

**Mean particle size:** 280 nm (250-320 nm)

**Hydration buffer:** 10% sucrose, 20mM HEPES, pH 7.5 ± 0.2

**Storage:** in refrigerator (2 - 8 degree C)

**Shelf-life:** 6 month for unopened vial

DOPC: 1,2-dioleoyl-sn-glycero-3-phosphocholine

CHOL: cholesterol

제조원 ; FormuMax Scientific Inc. (U.S.A)



삼보메디칼([www.sambomed.co.kr](http://www.sambomed.co.kr))

Tel: 02) 575-4945, Fax : 070) 4104-7595