

Placebo Control DOPC/CHOL/PEG-DSPE Liposomes

Code No. : F60203-P

Size : 1 ml & 5 ml

Product Description

Placebo Control DOPC/CHOL/PEG2000-DSPE Liposomes (100nm)

This placebo control liposomes is made to match the lipid composition and characteristics for the PEGylated fluorescent liposomes. The product is sterile filtered and filled in sterile vials guaranteed to be stable for 6 months when stored in refrigerator.

Specifications

Lipid composition: DOPC/CHOL/PEG2000-DSPE (55:45 mol/mol)

Lipid concentration: 50 mM (50-55 mM)

Dye: No dye

Mean particle size: 100 nm (85-110 nm)

Hydration buffer: 10% sucrose, 20mM HEPES, pH 7.5 ± 0.2

Stability: 6 months for unopened vials

Storage: 2-8 °C

DOPC: 1,2-dioleoyl-sn-glycero-3-phosphocholine

CHOL: cholesterol

제조원 ; FormuMax Scientific Inc. (U.S.A)



삼보메디칼(www.sambomed.co.kr)

Tel: 02) 575-4945, Fax : 070) 4104-7595