

Dilution Buffer for Doxoves™

Code No. : FB30020

Size : 2 ml & 10 ml

Image
Coming
Soon

Product Description

Dilution Buffer for Doxoves™

This dilution buffer made to match the bulk solution for Doxoves®

Specifications

Product code: FB30020

Composition: 10% sucrose, 10 mM histidine, pH 6.5 ± 0.2

Volume: 50 mL

Storage and Stability: sterile filtered (0.2um) , store at 2-8 degree C.

Shelf-life: 6 months stored in refrigerator. Use within 1 month weeks after first use.

Product For Research only and NOT for Human Use

제조원 ; FormuMax Scientific Inc. (U.S.A)



삼보메디칼(www.sambomed.co.kr)

Tel: 02) 575-4945, Fax : 02) 574-7595