

HRP conjugated rabbit anti-c-Myc

ORDERING INFORMATION

Catalog No: tAP-0002Size: $100 \mu l$ Storage: -20° C

Immunogen: c-Myc peptide
Ig Type: Rabbit IgG
Applications: WB & ELISA

Preparation: This antibody was produced by immunizing rabbit with KLH conjugated c-Myc peptide. The anti-c-Myc was affinity-purified using the peptide immobilized on a solid phase and conjugated with HRP.

Formulation and Storage: The purified HRP-conjugated anti-c-Myc was formulated in PBS containing 0.1% BSA, 50% glycerol and 0.01% thimerosal. The antibody should be kept at -20° C in a non auto-defrost freezer for long term storage.

Applications:

1. Western blot Yes (1:1,000-5,000)

2. ELISA Yes (1:5,000-20,000)

Optimal working dilutions for this anti-c-Myc should be determined for your particular assay conditions.

For Research Use Only

Contact & Ordering Information: Angio-Proteomie, 11 Park Drive, Suite 12, Boston, MA 02215, USA. Fax: (480) 247-4337, angioproteomie@gmail.com